

Fashion an Easier Lifestyle
with Assistive Technology



Personal Care

Long-Handled
Shoehorn



Sock Aid

Place aid between legs with cord towards floor.

Gather top of sock and line up sock heel with curve.



Sock Aid

Place aid on floor with curve side down.



Sock Aid

Pull cord up on both side of legs while inserting foot into sock.



Sock Aid

Continue to pull cords until sock is fully on foot and aid comes out of sock.



Furniture Riser



Plastic Food Guard



Suction Brush



Spoon and Fork Holders



Zipper and Button Puller



Food Preparation

Jelly Knife



Push/Pull Helper




Rocker Knife



Cutting Board

Suction cups hold board in place.



Jar Opener

Attach to wall. Screws included.



One-Hand Can Opener

Slide gray button back to release can opener.

Place opener on can.



One-Hand Can Opener

Squeeze opener until it snaps into place.

Continue squeezing and releasing opener until can is open.



One-Hand Can Opener

Slide back gray button to release can opener.



Housecleaning

Crate-On-Wheels



Long-Handled Toilet Brush



Sponge Wash Mitt



Recreation

Writing Guides Set



Hands-Free Card Holder




Standing Page Magnifier



Folding Book Stand



Pen/Pencil Cushions



Multi-Use

Reacher



Reacher Clips



Universal Cuff



Dycem Pads




Glow Tape

Combine with conventional lighting.
Remains visible 6-8 hours.



Foam



Seat Belt Handle



Gas Cap Wrench



Lamp Switch Enlargers



Key/Lock Identifiers



Thank You to Our Partners

Financial assistance provided by the Kansas Department of Health and Environment, Arthritis Program; the Arthritis Foundation, Kansas Chapter; and the Assistive Technology for Kansans Project.



This publication was supported by KDHE through USB/CCU722793-05 from Centers for Disease Control and Prevention (CDC). Its contents are solely the responsibility of the authors and do not necessarily represent the official views of CDC.



Kansas State University Agricultural Experiment Station and Cooperative Extension Service, Manhattan, KS.
Brand names appearing in this publication are for product identification purposes only. No endorsement is intended, nor is criticism implied of similar products not mentioned.

Contents of this publication may be freely reproduced for educational purposes. All other rights reserved. In each case, credit Debra M. Sellers and Kirsten Regan, *Fashion an Easier Lifestyle with Assistive Technology*, Kansas State University, July 2008.

For more information about adult development and aging, contact your local extension office. Kansas State University is an equal opportunity provider and employer.



www.aging.ksu.edu